


Paraffin-Embedding-Station



Modular Paraffin Embedding Station

Paraffin dispenser unit with 2 heatable trays for cassettes and / or molds, controllable paraffin outlet, Lighting, magnifying glass and cooling point.

Capacity Paraffin tank 4.5 liters.

Integrated weekly timer.

Technical specifications:


Dispenser unit: B560 x T640 x H400mm, 230V / 50Hz / 1KW

Max. 70°C,

Article number: AGS20



Big 4,5 Liter Paraffintank




Connection for a heatable forceps

Cassettes tray (160x190x40mm)

Illumination LED

Cooling Point

heated drawers



Optional:

Cooling unit for cooling the cassettes up to max. -13 ° C.

Cooling unit: B300 X T540 x H400mm, 230V / 50Hz / 0.4KW

Max. -13 ° C.

Article number: AGS20KP



Optional:


An additional heatable tray for accommodating the cassettes, large enough to accommodate most traditional tissue baskets.

Heating unit: B300 X T600 x H400mm, 230V / 50Hz / 0.4KW

Max. 70°C, capacity 2.5 liters

Article number: AGS20WP

Various heatable forceps



Cassettes Trim

Cassettes cleaning instrument



Cleans the cassettes of excess paraffin in order to enable an easy insertion into the clamp of the microtome.

Advantages:

No scratching with the knife more necessary.

No flying paraffin waste

No scratched numbers or barcodes

Technical specifications:

Weight: 2.7Kg

Dimensions: 240 x 160 x 170mm (D / W / H)

Temperature: max. 75 °C

230V 50Hz 0.1 KW

Article number: CT200