

Mobile Cooling Plates



Para Cooler TW **Mobile cooling plate with cooling tub**

Digital temperature controller.
Digital temperature display.
Wall and floor cooling.
Super Silent fans.
Temperature max.: -35°C .
264 paraffin blocks coolable at the same time.
Stainless steel cooling bath.

Technical specifications:

Dimensions of the device:
W / D / H: 485 x 610 x 670mm
Dimensions of the cooling surface:
W / D / H: 340 x 325 x 50mm
Refrigerant: 404A
Weight 42Kg
Connections: 230V 50Hz 0.4KW
Article No. 32002-7



Para Cooler TF **Mobile cooling plate flat**

Digital temperature display.
Super Silent fans.
Temperature max.: -18°C .
96 paraffin blocks can be cooled at the same time.
Cooling surface made of stainless steel.

Technical specifications:

Dimensions of the device:
W / D / H: 485 x 610 x 670mm
Dimensions of the cooling surface:
W / D / H: 345 x 344mm
Refrigerant: 134A
Weight 42Kg
Connections: 230V 50Hz 0.4KW
Article No. 32001-7



Table build in cooling plates

Para Cooler E

For installation in the table top
As a flat cooling plate or as a model with cooling tub.
Artikel Nr.: 30031 E
Artikel Nr.: 31000 E

Quality Made in Germany