

Slide and Cassette Printer

Cassette Printer SCP



Cassette printer for many standard cassettes.
The print is resistant to almost all common solvents.
Print from 2 lines up to 16 characters.
Automatic and semi-automatic counter.
Connection to a PC via USB port
Up to 10 freely programmable print images.

Article No. CP2008

Slide Printer



Slide Printer SSP

Object carrier printer for many standard commercial carriers.
The print is resistant to almost all common solvents.
Print 4 lines of up to 17 characters each.
Automatic and semi-automatic counter.
Connection to a PC via USB port.
Up to 10 freely programmable print images.

Technical specifications:
Dimensions: 370 x 220 x 150mm (W x D x H)
Power ratings: 230V 50Hz
Weight: 5Kg

Article No. SP2000

Printer ribbon for Vogel, Sakura und Surgipath printers.



Article No. 9000F