


MEPHISTO 12 Tissue Processor TP 12

MEPHISTO 14 Tissue Processor TP 14



Tissue processor with 12 or 14 stations.

Capacity: 100 or 300 cassettes.

Alternatively up to 1500 biopsies daily.

Up to 6 stations can be heated, except for paraffin also xylene and formalin are heated.

6 programs freely selectable and storable.




Agitation adjustable.

Weekend program with up to 99 hours start delay.

Emergency power supply to protect the Tissue samples.

Can be connected to a laboratory exhaust.

No valves, no clogged pipes, no cleaning fluids, no cleaning program, no maintenance costs, intervention in the current program is possible at any time.



Intuitive user interface.



Easy solvent change.
The front panel can be opened completely.

Technical specifications:
Dimensions: B1100 x D 630 x H 769mm
Dimensions: B 895 x D 630 x H 769mm
Reagent container: stainless steel 3.5 liters
Power supply: 230V / 50Hz / 1,2KW
Weight: 66 Kg / 77Kg

Quality Made in Germany