

Fast Fixation devices.

Fast fixation instrument Quick-Fix M

Fast fixation instrument for big tissue samples.

Fast fixation of tissue samples with a diameter of up to 20cm.
Excellent fixation in a short time.
Without microwaves.
Temperature from room temperature to 60 °C free selectable.
Automatic pressure / vacuum system.
Continuous centrifugal rinsing of the tissue.
Small particulate filter.
Retort with a capacity of 15 L of formalin.
Several tissue samples can be fixed at the same time.

Special basket to hold large
Tissue samples, thickness of 5mm - 10mm.
Technische Daten

Dimensions (W / L / H): 500x420x520 mm
Weight: 30 kg
Power ratings: 230 V; 50 Hz; 0.8 KW (CE)
Article No .: 40070



Fast fixation instrument Quick-Fix

Fast fixation instrument for small and medium tissue samples.

Fast fixation device with integrated tissue filter and automatic suction of the formalin vapors
On request with downstream activated carbon filter.
The fixing time is reduced to a minimum.
Excellent for fat tissue.
Digital temperature control.
Capacity: 75 cassettes.

Technical data:
Dimensions of the device (L / W / H) 435 x 290 x 440 mm
Weight: 22 Kg
Power ratings: 230V, 50Hz, 0.3 KW
Article No .: 40050



Qualität Made in Germany