

Fast freezing unit Cryofix -45 °C



Instrument for fast freezing of tissue samples.

Fits all conventional cryostats.

Absolute alcohol is present in the retort up to -45°C
The formation of ice crystals almost becomes completely avoided.
Immediate shock freezing reduces the freezing time to a minimum.
No danger as with liquid nitrogen.

The compact design allows installation close to your Cryostat.
Also as an alternative to the formalin fixation in sensitive Antigens.

Technical specifications:
WxDxH: 300x410x460mm
Weight: 31Kg
Power supply: 230V / 50Hz / 0.3KW
Article No. 038001

Incubator / slides dryer

Our smallest fits in every laboratory



Temperature: room temperature to + 70 °C
Temperature setting: +/- 1 degree
Built-in fan
Actual value Display on the LCD display
Stainless steel inner body
2 stainless steel shelves removable / height adjustable
Stainless steel paraffin collecting basin

Technical specifications:
Overall dimensions: 340x345x525 mm (WxDxH)
Internal dimensions: 230x275x325 mm (WxDxH)
Operating voltage: 230 V / 50Hz
Power consumption: 100 VA
Weight: 28 kg
Article 90002

Qualität Made in Germany