

Sterilizer / autoclave for immunohistology

Controlled demasking of paraffin sections.

For reproducible results



The electronic autoclave with stainless steel chamber

Electronic autoclave from room temperature to 140°C.
Several programs preinstalled.

Controller unit has a printer port for the autoclave receipt printer and a USB port for PC connection.

Basic version of PC software for data monitoring included (professional license required for advanced programming functions).

On request fractional flow method.

Technical specifications:

12 liter device

24cm
26cm
35cm

diameter inside
height loadable
total height

18 liter device

24cm
39cm
46.5cm

Power ratings: 230V / 50Hz / 1.9KW

Filling



Locking



Start



Qualität Made in the EU