

Waterbath with illumination



Round water bath with illumination
Digital temperature display and control.
Water bath removable.
Only the glass pan is transported for water change.
Cleaning of the glass pan in the dishwasher.
Side shelf for slides.
Water bath: Diameter 210mm height 45mm
Outside dimensions: W310 x L355 x H170mm
Temperature max. 60 ° C
Power ratings: 230V / 50Hz / 200W
Article No. WBRL



Square water bath with illumination
Digital temperature display and control.
Water bath removable.
Only the glass pan is transported for water change.
Cleaning the glass pan in the dishwasher.
Side shelf for slides.
Water bath: 210 x 210mm Height 55mm
Overall dimensions: W310 x L355 x H170mm
Temperature max. 60 ° C
Power ratings: 230V / 50Hz / 200W
Article No. WBQL

Determine the contrast by selecting the appropriate color.
The lighting is dimmable and allows fatigue-free working.



Water bath heated plate combination with lighting



Square or round water bath heated plate with illumination
Digital temperature display and control.
Water bath removable

Technical specifications:

Water bath: 210 x 210 x 55mm
Overall dimensions: 470 x 295 x 94mm
Heated plate: 270 x 180mm
Weight: 6.7 Kg
Capacity: 2.0 L
Temperature max. 60 ° C
Power ratings: 230V / 50Hz / 400W
Article No. WBST 2000

Quality Made in Germany