

## Waterbath round



Round water bath with extra wide edge.  
Splash-proof on / off switch.  
Heating control light.

Technical specifications:  
Internal dimensions: 200 x 60mm  
Outside dimensions: 330 x 110mm  
Weight: 2.7 kg  
Capacity: 1.7L  
Temperature max. 90 °C  
Power ratings: 230V / 50Hz / 250W  
Article No. WBR

## Heated Plate



Slide drayer  
Stretching bench

Digital temperature display.  
Digital temperature control.  
Working surface only 45mm high.

Technical specifications:  
Dimensions: 550 x 180 x 45 (L x W x H)  
Working area: 155 x 450 mm.  
Working height: 45 mm  
Temperature range: 30°C - 80°C  
Power ratings: 230V / 50Hz / 300W.

**Quality Made in Germany**

## Paraffin Dispenser



Capacity: 15 liters  
On request with built in filter for melting  
of cutting waste  
Temperature of 20 °C - 70 °C continuously adjustable.  
Outlet faucet heated.  
PTC heater with overheat protection.  
Stainless steel lid.

### Technical specifications:

Capacity: max. 15L  
B / T / H: 380x420x520  
Weight: 14 Kg  
Power ratings: 230V / 50Hz / 550W  
Article No. 70010

**Quality Made in Germany**