



**Ralf W. Weinkauff**  
Medizin und Umwelttechnik

Ahornweg 4a  
91301 Forchheim

Tel. 0 91 91/970 661

Fax 0 91 91/970 665

[www.weinkauff-medizintechnik.de](http://www.weinkauff-medizintechnik.de)

email: [info@weinkauff-medizintechnik.de](mailto:info@weinkauff-medizintechnik.de)



## Rotary Microtom CUT 4060 E

Quality Made in Germany



Hand Control



Optional  
Foot Pedal

Range of section thickness setting from 0,5  $\mu\text{m}$  to 60  $\mu\text{m}$

Section thickness feed 0,5  $\mu\text{m}$  from 0,5  $\mu\text{m}$  to 2  $\mu\text{m}$

1  $\mu\text{m}$  from 2  $\mu\text{m}$  to 10  $\mu\text{m}$

2  $\mu\text{m}$  from 10  $\mu\text{m}$  to 20  $\mu\text{m}$

5  $\mu\text{m}$  from 20  $\mu\text{m}$  to 60  $\mu\text{m}$

Horizontal Specimen Advance 40 mm

Vertikal Stroke 70 mm

Specimen Size Standard Maximal 55 mm x 55 mm

Dimension Length x Width x Height:

550 mm x 460 mm x 280 mm

Weight (without accessories) 45 kg

Manufactured according to the standards of class 4/MedGV.

The best choice for research and routine work in biological, medical or industrial areas.

Guaranteed reliability through over-extension protection and sliding-clutch anti-locking system in the coarse advance.

Highest precision of sections from paraffin or plastic blocks over the total advance range.