



**Ralf W. Weinkauff**  
Medizin und Umwelttechnik

Ahornweg 4a  
91301 Forchheim

Tel. 0 91 91/970 661  
Fax 0 91 91/970 665

www.weinkauff-medizintechnik.de  
email: info@weinkauff-medizintechnik.de



**Rotary Microtom CUT 4060**



**Rotary Microtom CUT 4060 F**

**Rotary Microtom CUT 4060 Technical Details:**

Section Thickness range: 05 - 60 $\mu$ m

Thickness increments: 0 - 60 $\mu$ m in 0,5 increments

Maximum section area: 55 x 50mm

Maximum horizontal specimen tilt: 40mm

Maximum vertical specimen tilt: 70mm

Movable knifeholder

Security handwheel lock

Included:

Knifeholder for standart knife or disposable blades.

Standart bloc holder or cassettes holder.

Weight: 32 Kg

Made in Germany