

Slides & Paraffin Block Archives

Drawer cabinets with central lock



A. Storage cabinet for storage of approx. 62,000 standard specimen slides 76 x 26 mm or about 16,800 paraffin blocks with double storage
6 full extension drawers (100% extendable) incl. Trays
Dimensions: H x W x D 850 x 745 x 750 mm




B. Storage cabinet for the storage of approx. 92,880 standard slide carriers 76 x 26 mm or approximately 25,200 paraffin blocks with double storage
9 full extension drawers (100% extendable) incl. Trays
Dimensions: H x W x D 1225 x 745 x 750 mm



C. Storage cabinet for storage of approx. 148,000 standard slide carriers 76 x 26 mm or approximately 33,600 paraffin blocks with double storage
12 full extension drawers (100% extendable) incl. Trays
Measurements: H x W x D 1600 x 745 x 750 mm

D. Storage cabinet for storage of approx. 92,880 standard slide carriers 76 x 26 mm or approximately 25,200 paraffin blocks with double storage
6 full extension drawers (100% extendable) incl. Trays
Dimensions: H x W x D 850 x 1046 x 750 mm

Slide Trays



F. Storage cabinet for the storage of approx. 139.320 standard slide carriers 76 x 26 mm or approximately 37,800 paraffin blocks with double storage
6 full extension drawers (100% extendable) incl. Trays
Dimensions: H x W x D 1225 x 1046 x 750 mm

G. Storage cabinet for storage of approx. 185.760 standard slide carriers 76 x 26 mm or approximately 50,400 paraffin blocks with double storage
6 full extension drawers (100% extendable) incl. Trays
Measurements: H x W x D 1600 x 1046 x 750 mm

Paraffinblock Tray double



I. Storage cabinet for the storage of approx. 154,000 standard slide carriers 76 x 26 mm or about 21,200 paraffin blocks with simple storage
10 full extension drawers (100% extendable) incl. Trays
Measurements: H x W x D 1100 x 1046 x 750 mm

Safety:

All cabinets are lockable.

No risk of tipping, only one drawer can be opened at any time.

Advantages:

Optimum use of space, paraffin blocks are stored in trays double-storey.

It is possible to store mixed slides and paraffin blocks in a cabinet.

Delivery to the place of use and assembly by our professional staff possible.

Qualität Made in Germany