

Dimensions of the appliance:
W/D/H: 415 x 460 x 300mm
Dimensions of the cooling surface:
W/D 345 x 344mm

Table Top Cooling Plate

Blue Cooler the energy-saving cooling plate.

The largest cooled surface with the least space requirement.
Digital temperature regulation
Super Silent fans.
Temperature max.: -18°C .
96 paraffin blocks can be cooled at the same time.
Stainless steel frame.
Cooling surface made of stainless steel.


Technical specification:
Refrigerant: 600A
Weight 14Kg
Power supply 230V/50Hz/300W
Article No. 31400



Para Cooler O Table cooling plate with cooling tub.

Digital temperature regulation
Wall and floor cooling.
Super Silent fans.
Temperature max.: -35°C .
264 paraffin blocks coolable at the same time.
Stainless steel frame.
Stainless steel cooling bath.

Technical specifications:
Dimensions of the device:
W / D / H: 415 x 600 x 320mm
Dimensions of the cooling surface:
W / D / H: 340 x 325 x 50mm
Refrigerant: R 290
Weight 24Kg
Power ratings: 230V / 50Hz / 300W
Article No. 30031



Para Cooler A cooling plate flat

Digital temperature regulation
Super Silent fans.
Temperature max.: -18°C .
96 paraffin blocks can be cooled at the same time.
Stainless steel frame.
Cooling surface made of stainless steel.

Technical specifications:
B / T / H: 415 x 600 x 280mm
Dimensions of the cooling surface:
B / T 345 x 344mm
Refrigerant: 600A
Weight 18Kg
Power ratings: 230V / 50Hz / 300W
Article 31000

Quality Made in Germany