


Linear Stainer LS 72 for standard und special stains.

Multi stainer with 72 cuvettes

Multiple staining with up to 40 slide baskets simultaneously possible.



Staining machine for highest demands with a maximum capacity of 1,800 slides per hour.

Up to 20 programs freely selectable.

Loading station for 4 slide baskets.

Unloading station for 4 slide baskets.

With up to 6 fresh water cuvettes.

Agitation adjustable.

Optional built-in activated carbon filter or connection to an existing exhaust air system.

Adaptable to almost all automatic coverslippers.

Slide holder for 20 or 30 slides.

Technical data:

Dimensions: W x D x H: 200 x 75 x 120cm

Connected load: 230V / 50Hz / 1,2KW

Weight: 210 Kg Article No. 72000