


Table cooling plates



Technical specifications:
Dimensions of the device:
W / D / H: 415 x 460 x 300mm
Dimensions of the cooling surface:
W / D 345 x 344mm

Blue Cooler the energy-saving cooling plate

The largest cooling surface with the smallest footprint.

Energy consumption only 0.07 KW at -15 ° C

Digital temperature control

Super silent fan.

Temperature max.: -18 ° C.

96 paraffin blocks can be cooled simultaneously.

Frame made of stainless steel.

Cooling surface made of stainless steel.

Refrigerant: R600a

Weight 14Kg

Connected load: 230V / 50Hz / 300W

Article No. 31400



Para Cooler O Tischkühlplatte mit Kühlwanne.

Para Cooler O table cooling plate with cooling tub.

Digital temperature control

Wall and floor cooling.

Super silent fan.

Temperature max.: -35 ° C.

264 paraffin blocks can be cooled simultaneously.

Frame made of stainless steel.

Cooling tub made of stainless steel.

Technical specifications:

Dimensions of the device:

W / D / H: 415 x 600 x 320mm

Dimensions of the cooling surface:


W / D / H: 340 x 325 x 50mm

Refrigerant: R290

Weight 32Kg

Connected load: 230V / 50Hz / 300W

Article No. 30031



Para Cooler A flat table cooling plate.

Digital temperature control

Super silent fan.

Temperature max .: -18 ° C.

96 paraffin blocks can be cooled simultaneously.

Frame made of stainless steel.

Cooling surface made of stainless steel.

Technical specifications:

Dimensions of the device:

W / D / H: 415 x 600 x 280mm

Dimensions of the cooling surface:

B / D 345 x 344mm

Refrigerant: R600a

Weight 24Kg

Connected load: 230V / 50Hz / 300W

Article No. 31000

Quality Made in Germany