



Electrically heated tweezers.



Suitable for our tweezers station.
It can connect tweezers in sizes from 0.5mm to 4mm.
The individual tweezers are color-coded according to size for easy orientation, marked. All tweezers are heated only at the tip to approx. 68 degrees.
at the tip only. The connecting cable is heat resistant up to 80 degrees, highly flexible, and resistant to kerosene.

Order numbers:



4709205 (0,5 mm tweezers) blue

4709210 (1 mm tweezers) red

4709220 (2 mm tweezers) yellow

4709240 (4 mm tweezers) green

All tweezers can also be connected to our embedding center

Electrical connection: 24V DC 0,34A