

## Formalin filling device

Filling device for small or large containers



The liquid dispenser with is the optimal solution for filling formalin in containers in any pathology lab. Due to the multiple adjustment options the desired amount can be filled in containers, from 1ml up to 1000ml is possible

Advantages:

- Integrated exhaust system with activated carbon filter.
- Filtering by activated carbon filter plates, no granules, thereby Reduction of fine dust pollution.
- Splash-proof dosing switch
- Electronic filling control from 1ml - 1000ml.
- Stainless steel check valve with built-in filter



On request with  
a filling gun.

Technical specifications:

Dimensions of the device (L / W / H) 300 x 220 x 160mm

Weight: 6,5 Kg

Power requirements: 230V, 50Hz, 0,2 KW

Article Nr. 50010