

Histology Wax Recycling

Paraffin Recycler

Reprocessing of with xylene and alcohol contaminated Paraffin from Tissue Processors.

Recycling of old Paraffin blocks.

15 liter capacity.

Sophisticated vacuum system down to minus 0.8 bar.

Saving of storage and disposal costs.

Integrated overheating protection.



Almost 100% paraffin recovery through triple filtration.



Technical specifications:

Dimensions of the device: L/W/H 500x420x520mm

Weight 30kg

230V/50Hz/800W

Article No. 60010