

# Waterbath round



## Design and Safety Features

The round water bath is engineered with an extra wide edge, which provides increased stability during operation and makes the unit easier to handle. This thoughtful design ensures users can confidently transport and position the water bath as needed.

For enhanced safety, the water bath features a splash-proof on/off switch. This important addition helps protect the electrical components from moisture, reducing the risk of electrical faults or accidents during use.

In addition, the water bath is equipped with a heating control light. This indicator allows users to easily observe whether the heating element is currently active, offering a clear visual cue regarding the operational status of the unit.

## Technical Data

Internal dimensions	Outside dimensions	Weight	Capacity	Temperature max.	Power ratings	Article No.
200 x 60mm	330 x 110mm	2.7 kg	1.7 L	90°C	230V / 50Hz / 250W	WBR