


# Waterbath with illumination



## Features

The water bath is designed with a practical square shape, allowing for efficient use of space within the laboratory environment. Integrated illumination enhances visibility, making it easier to observe the contents during operation.

User-friendly digital controls are provided, enabling precise adjustment of settings. A temperature display is included to ensure accurate and convenient monitoring of the water bath's operating conditions.

For added convenience, the unit features designated storage areas for slides located both at the front and the rear. This thoughtful design detail supports efficient workflow and organization during laboratory procedures.

## Technical Data

<b>Water bath dimensions</b>	<b>260 x 250mm height 60mm</b>
<b>Outside dimensions</b>	B320 x L440 x H90mm
<b>Temperature max.</b>	60 °C
<b>Capacity</b>	2.5 liter
<b>Weight</b>	5.1 Kg
<b>Connected load</b>	230V / 50Hz / 500W
<b>Article No.</b>	WBQL2.5