

# Liquid Dispenser

for non flammable liquids



## Key Features of the Liquid Dispenser

The liquid dispenser is equipped with an integrated exhaust system, making it one of the leading options for bottling non-flammable liquids in modern laboratory environments. The dispenser allows for precise volume adjustment, enabling users to fill bottles or liquid containers of various sizes efficiently. Additionally, it supports the pre-filling of formalin containers, enhancing workflow efficiency.

## Operator Safety

To ensure user safety, the dispenser includes features designed to protect operators from exposure to harmful vapors. This is achieved through an exhaust system that incorporates an activated charcoal filter, which effectively captures and neutralizes hazardous emissions.

## Precision and Control

The dispenser is equipped with electronic bottling controls that provide adjustable dispensing volumes ranging from 1 ml to 1000 ml. This allows for accurate and customizable filling, accommodating a wide spectrum of laboratory requirements.

## Durability and Optional Accessories

A stainless steel check valve with an integrated filter is included to ensure reliable operation and maintain the integrity of dispensed liquids. For added convenience, an optional hand drainage faucet is also available.



Dimensions (L/W/H)	Weight	Power supply
300 x 220 x 160 mm	3,5 Kg	230V, 50Hz, 0,2 KW