

Electric Heated Tweezers



Order No. 4509200S2

Compatibility and Features

These electric heated tweezers are designed to be fully compatible with our tweezers station. They support a range of tweezer sizes, accommodating tips from 0.5 mm to 4 mm. Each tweezer is color-coded by size, ensuring effortless identification and orientation during use.

Heating and Safety Specifications

The heating element is concentrated at the tip of each tweezer, maintaining a consistent operating temperature of approximately 68 degrees. The connecting cable is engineered for durability, being both highly flexible and heat resistant up to 80 degrees, which ensures safe and reliable operation.

Product Variants

Order Number	Size	Color Code	
4709205	0.5 mm	Blue	
4709210	1 mm	Red	
4709220	2 mm	Yellow	
4709240	4 mm	Green	

Additional Compatibility

All tweezers are also suitable for connection to our embedding center, expanding their versatility within laboratory workflows.

Electrical Specifications

The tweezers operate on a 24V DC power supply with a current of 0.34A.